### Produktbeschreibung

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### Product description

CATALYST UNI BOND is a product for additional and/or extended self etch applications of universal dental adhesive UNI BOND for all the cases, where light cure of the adhesive is difficult or impossible (e.g. in the root canal).

The mixture of UNI BOND/ CATALYST UNI BOND is suitable for moist, wet and dry surfaces.

### Indications/Intended use

Activator for the self-etch technique for:

- adhesive cementing of posts and subsequent core build-up with dual or self cure composites.

### Performance features

The performance features of the product meet the requirements of the intended use.

### Contraindications

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material.

### Patient target group

Persons who are treated during a dental procedure.

### Intended users

This medical device should only be used by a professionally trained dental practitioner.

### Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the material. Neither store the material in proximity of eugenol containing products, nor let the material allow coming into contact with materials containing eugenol.

### Application

#### 1. Isolation

Rubber dam is the recommended method of isolation.

#### 2. Root Canal Preparation

Prepare the root canal appropriate to the selected post.

#### 3. Application of the Adhesive

One drop of UNI BOND and one drop of CATALYST UNI BOND were combined in a mixing pallet and mixed for **5-10 seconds** under subdued light.

#### Notes:

Do not interchange lids of the bottles, because this can lead to a cross-contamination of the liquids.

Apply the UNI BOND/CATALYST UNI BOND mixture with a suitable brush onto the root canal walls for 30 seconds with agitation. Applying the adhesive on the surface without agitation is inadequate. Repeat procedure 1–2 times. The material should build a homogeneous layer. Remove excess material thoroughly e.g. with paper points.

Air thin gently (**at least 5-10 seconds**) to remove the volatile components and to disperse the adhesive until a glossy and uniform layer results. For optimal adhesion light cure with a dental halogen light unit or an LED (wavelength 400–500 nm, light intensity min. 1000 mW/cm<sup>2</sup>) for **20 seconds** is recommend.

#### 4. Post Cementing

Prepare the post according to manufacturer instructions.

Place the post with a flowable dual cure or self cure composite into the root canal.

Refer for this and for core build-up to manufacturer instructions of the dual or self cure composite.

#### Additional informations/Warnings

- Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates.
- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

#### Storage

Do not store above 25 °C (77 °F)! Avoid storage in direct sunlight. Do not use after expiration date.

#### Disposal

Disposal of the product according to local authority regulations.

#### Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

#### Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior Dental Materials GmbH product.

#### Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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P.L. Superior Dental Materials GmbH, Stiller Weg 15 A, 22607 Hamburg, Germany